

## Surgical face mask Type IIR (BFE ≥ 98%) Kids



### From 6-9 years

#### Información técnica

Typology: Disposable face mask for kids

Model: 150 mm x 65 mm

Standard: EN 14683-2019 + AC:2019. Type IIR

Size: 150 mm x 65 mm (6-9 years) Bacterial filtration efficiency: ≥ 98%

Features:

Type IIR: maximum surgical quality


Splash resistant

Structure: Mask body, nose clip and mask belt

Validity period: 24 months

Allergy warning: Non-woven textile

### Non-woven spunbond polypropylene fabric



Melt-blown layer with 99% BFE

### **Fitting instructions**

- 1. With the blue side facing outwards, place it over the nose and mouth, with the nose clip facing upwards. Then place the loops behind the ears.
- 2. Gently press on the nose clip in order to adjust it to the nose.
- 3. Pull the lower part of the mask so that it is beneath the chin and check that it fits well.







#### Important information

- Hypo-allergenic waterproof polypropylene material, including flexible nasal clip.
- CE Test Report issued by TÜV SUD.
- Conditions for storage: To be stored in an area with relatively low humidity, free from corrosive gases and well ventilated, avoiding high temperatures (>40°C).



# Surgical face mask Type IIR (BFE ≥ 98%) Kids

### From 6-9 years









# Surgical face mask Type IIR (BFE ≥ 98%) Kids

From 6-9 years



CHARACTERISTICS	STANDARD	RESULTS
Bacterial Filtration Efficiency (BFE)	EN 14683:2019+AC:2019 (E) Annexe B	BFE ≥ 99%
Differential pressure (ΔP)	EN 14683:2019+AC:2019 (E) Annexe C	ΔP < 38 Pa/cm2
Microbial cleaning	EN 14683:2019+AC:2019 (E) Annexe D	CFU/g < 20
Synthetic blood penetration	ISO 22609:2004	Not detected

Source: Test Report No.: 200193327 (TÜV SÜD)

### 10 pcs bag specifications

Quantity: 10 units/bag

Dimensions: 210 x 120 x 23 mm

Weight: 25 g

SKU: MASK-ST-00025.0010

6 - 9 years



### **Master carton specifications**

Quantity: 100 bags/master carton

Total quantity: 1000 units

Dimensions: 520 x 370 x 250 mm

Weight: 4.8 kg

